

BFAR Rabbit Polyclonal Antibody

Catalog No: #55475

Package Size: #55475-1 50ul #55475-2 100ul


Description

Product Name	BFAR Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human BFAR (NP_057645.1).
Conjugates	Unconjugated
Other Names	BFAR;BAR;RNF47
Accession No.	Swiss Prot:Q9NZS9GeneID:51283
Calculated MW	38kDa/52kDa
SDS-PAGE MW	53kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100


Images




Western blot analysis of extracts of various cell lines, using BFAR at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using BFAR Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH-3T3 cells using BFAR Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of U-2 OS cells using BFAR Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.