

NID1 Rabbit Polyclonal Antibody

Catalog No: #55485

Package Size: #55485-1 50ul #55485-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NID1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human NID1 (NP_002499.2).
Conjugates	Unconjugated
Other Names	NID1;NID;nidogen-1
Accession No.	Uniprot:P14543GenelD:4811
Calculated MW	122kDa/136kDa
SDS-PAGE MW	178kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NID1 antibody.

Background

This gene encodes a member of the nidogen family of basement membrane glycoproteins. The protein interacts with several other components of basement membranes, and may play a role in cell interactions with the extracellular matrix.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.