

PGLS Rabbit Polyclonal Antibody

Catalog No: #55493

Package Size: #55493-1 50ul #55493-2 100ul


Description

Product Name	PGLS Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PGLS (NP_036220.1).
Conjugates	Unconjugated
Other Names	PGLS;6PGL;HEL-S-304
Accession No.	Swiss Prot:O95336GeneID:25796
Calculated MW	27kDa
SDS-PAGE MW	28kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PGLS at 1:1000 dilution.

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.