

RAB24 Rabbit Polyclonal Antibody

Catalog No: #55497

Package Size: #55497-1 50ul #55497-2 100ul


Description

Product Name	RAB24 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RAB24 (NP_001026847.1).
Conjugates	Unconjugated
Other Names	RAB24
Accession No.	Uniprot:Q969Q5GeneID:53917
Calculated MW	23kDa
SDS-PAGE MW	23kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RAB24 antibody.

Background

RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Oikkonen et al., 1993 [PubMed 8126105]).[supplied by OMIM, Aug 2009]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.