

TGFB3 Rabbit Polyclonal Antibody

Catalog No: #55521

Package Size: #55521-1 50ul #55521-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TGFB3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TGFB3 (NP_003230.1).
Conjugates	Unconjugated
Other Names	TGFB3;ARVD;ARVD1;LDS5;RNHF;TGF-beta3
Accession No.	Uniprot:P10600GenelD:7043
Calculated MW	35kDa/47kDa
SDS-PAGE MW	60kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TGFB3 antibody.

Background

This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and is involved in embryogenesis and cell differentiation. Defects in this gene are a cause of familial arrhythmogenic right ventricular dysplasia 1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.