

# N4BP1 Rabbit Polyclonal Antibody

Catalog No: #55525

Package Size: #55525-1 50ul #55525-2 100ul



**SAB**<sup>®</sup>  
Signalway Antibody

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                          |
|-----------------------|----------------------------------------------------------|
| Product Name          | N4BP1 Rabbit Polyclonal Antibody                         |
| Host Species          | Rabbit                                                   |
| Clonality             | Polyclonal                                               |
| Isotype               | IgG                                                      |
| Purification          | Affinity purification                                    |
| Applications          | WB,IF                                                    |
| Species Reactivity    | Human,Mouse                                              |
| Immunogen Description | Recombinant fusion protein of mouse N4BP1 (NP_085040.2). |
| Conjugates            | Unconjugated                                             |
| Other Names           | N4BP1                                                    |
| Accession No.         | Swiss Prot:O75113GeneID:80750                            |
| Calculated MW         | 99kDa                                                    |
| SDS-PAGE MW           | 120kDa                                                   |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.  |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.              |


## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using N4BP1 at 1:1000 dilution.



Immunofluorescence analysis of U2OS cells using N4BP1 at dilution of 1:100. Blue: DAPI for nuclear staining.

## Background

---

–

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.