

FCGRT Rabbit Polyclonal Antibody

Catalog No: #55536

Package Size: #55536-1 50ul #55536-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FCGRT Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human FCGRT (NP_004098.1).
Conjugates	Unconjugated
Other Names	FCGRT;FCRN;alpha-chain
Accession No.	Swiss Prot:P55899GenelD:2217
Calculated MW	39kDa
SDS-PAGE MW	55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using FCGRT at 1:1000 dilution.

Background

This gene encodes a receptor that binds the Fc region of monomeric immunoglobulin G. The encoded protein transfers immunoglobulin G antibodies from mother to fetus across the placenta. This protein also binds immunoglobulin G to protect the antibody from degradation. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.