

AIFM3 Rabbit Polyclonal Antibody

Catalog No: #55554

Package Size: #55554-1 50ul #55554-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AIFM3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human AIFM3 (NP_001018070.1).
Conjugates	Unconjugated
Other Names	AIFM3;AIFL
Accession No.	Swiss Prot:Q96NN9GeneID:150209
Calculated MW	65kDa/66kDa
SDS-PAGE MW	72kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat brain using AIFM3 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded rat kidney using AIFM3 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human lung cancer using AIFM3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver cancer using AIFM3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using AIFM3 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using AIFM3 at 1:1000 dilution.

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.