

URG11 Rabbit Polyclonal Antibody

Catalog No: #55555

Package Size: #55555-1 50ul #55555-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	URG11 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human URG11 (NP_689931.2).
Conjugates	Unconjugated
Other Names	VWCE;URG11;VWC1
Accession No.	Swiss Prot:Q96DN2Gene ID:220001
Calculated MW	22kDa/99kDa
SDS-PAGE MW	100kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB□1:2000 - 1:5000IHC□1:50 - 1:200

Images



Western blot analysis of extracts of HepG2 cells, using VWCE Antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.