

# GMFB Rabbit Polyclonal Antibody

Catalog No: #55563

Package Size: #55563-1 50ul #55563-2 100ul


## Description

|                       |                                                         |
|-----------------------|---------------------------------------------------------|
| Product Name          | GMFB Rabbit Polyclonal Antibody                         |
| Host Species          | Rabbit                                                  |
| Clonality             | Polyclonal                                              |
| Isotype               | IgG                                                     |
| Purification          | Affinity purification                                   |
| Applications          | WB,IF                                                   |
| Species Reactivity    | Human,Mouse,Rat                                         |
| Immunogen Description | Recombinant fusion protein of human GMFB (NP_004115.1). |
| Conjugates            | Unconjugated                                            |
| Other Names           | GMFB;GMF                                                |
| Accession No.         | Swiss Prot:P60983GenelD:2764                            |
| Calculated MW         | 16kDa                                                   |
| SDS-PAGE MW           | 17kDa                                                   |
| Formulation           | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage               | Store at -20°C. Avoid freeze / thaw cycles.             |

## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using GMFB at 1:1000 dilution.

## Background

—

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.