

PFDN1 Rabbit Polyclonal Antibody

Catalog No: #55564

Package Size: #55564-1 50ul #55564-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PFDN1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PFDN1 (NP_002613.2).
Conjugates	Unconjugated
Other Names	PFDN1;PDF;PFD1
Accession No.	Swiss Prot:O60925GeneID:5201
Calculated MW	14kDa
SDS-PAGE MW	14kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using PFDN1 at 1:3000 dilution.

Background

This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.