

ANXA3 Rabbit Polyclonal Antibody

Catalog No: #55568

Package Size: #55568-1 50ul #55568-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ANXA3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ANXA3 (NP_005130.1).
Conjugates	Unconjugated
Other Names	ANXA3;ANX3
Accession No.	Uniprot:P12429GenelD:306
Calculated MW	36kDa
SDS-PAGE MW	36kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ANXA3 antibody.

Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions in the inhibition of phospholipase A2 and cleavage of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. This protein may also play a role in anti-coagulation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.