

SCFD1 Rabbit Polyclonal Antibody

Catalog No: #55574

Package Size: #55574-1 50ul #55574-2 100ul


Description

Product Name	SCFD1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SCFD1 (NP_057190.2).
Conjugates	Unconjugated
Other Names	SCFD1;C14orf163;RA410;SLY1;SLY1P;STXBP1L2
Accession No.	Swiss Prot:Q8WVM8GeneID:23256
Calculated MW	62kDa/65kDa/72kDa
SDS-PAGE MW	72kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using SCFD1 at 1:1000 dilution.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.