

GPM6A Rabbit Polyclonal Antibody

Catalog No: #55585

Package Size: #55585-1 50ul #55585-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GPM6A Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human GPM6A (NP_005268.1).
Conjugates	Unconjugated
Other Names	GPM6A;GPM6;M6A
Accession No.	Swiss Prot:P51674GenelD:2823
Calculated MW	29kDa/30kDa/31kDa
SDS-PAGE MW	31kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of U-87MG cells, using GPM6A at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded mouse stomach using GPM6A at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using GPM6A at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.