

SH3PXD2A Rabbit Polyclonal Antibody

Catalog No: #55596

Package Size: #55596-1 50ul #55596-2 100ul


Description

Product Name	SH3PXD2A Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human SH3PXD2A
Conjugates	Unconjugated
Other Names	SH3PXD2A;FISH;SH3MD1;TKS5
Accession No.	Swiss Prot:Q5TCZ1Gene ID:9644
Calculated MW	102kDa/121kDa/125kDa
SDS-PAGE MW	100kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of HeLa cells, using SH3PXD2A antibody.



Immunofluorescence analysis of L929 cells using SH3PXD2A Rabbit pAb.

Background

—

Note: This product is for in vitro research use only and is not intended for use in humans or animals.