

NUP50 Rabbit Polyclonal Antibody

Catalog No: #55597

Package Size: #55597-1 50ul #55597-2 100ul


Description

Product Name	NUP50 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NUP50 (NP_009103.2).
Conjugates	Unconjugated
Other Names	NUP50;NPAP60;NPAP60L
Accession No.	Swiss Prot:Q9UKX7GeneID:10762
Calculated MW	46kDa/50kDa
SDS-PAGE MW	58kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000 IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using NUP50 at 1:3000 dilution.

Background

The nuclear pore complex is a massive structure that extends across the nuclear envelope, forming a gateway that regulates the flow of macromolecules between the nucleus and the cytoplasm. Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells. The protein encoded by this gene is a member of the FG-repeat containing nucleoporins that functions as a soluble cofactor in importin- α : β -mediated nuclear protein import. Pseudogenes of this gene are found on chromosomes 5, 6, and 14. Two transcript variants

encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.