

SLC30A1 Rabbit Polyclonal Antibody

Catalog No: #55639

Package Size: #55639-1 50ul #55639-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC30A1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human SLC30A1 (NP_067017.2).
Conjugates	Unconjugated
Other Names	SLC30A1;ZNT1;ZRC1
Accession No.	Uniprot:Q9Y6M5GeneID:7779
Calculated MW	55kDa
SDS-PAGE MW	55kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse kidney, using SLC30A1 antibody.

Background

May be involved in zinc transport out of the cell.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.