

RRAGD Rabbit Polyclonal Antibody

Catalog No: #55644

Package Size: #55644-1 50ul #55644-2 100ul


Description

Product Name	RRAGD Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RRAGD (NP_067067.1).
Conjugates	Unconjugated
Other Names	RRAGD;RAGD;bA11D8.2.1
Accession No.	Swiss Prot:Q9NQL2GeneID:58528
Calculated MW	28kDa/45kDa
SDS-PAGE MW	46kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:1000 - 1:4000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using RRAGD at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded human liver cancer using RRAGD at dilution of 1:100 (40x lens).

Background

-

Note: This product is for in vitro research use only and is not intended for use in humans or animals.