

Biotin Rabbit pAb

Catalog No: #55648

Package Size: #55648-1 50ul #55648-2 100ul


Description

Product Name	Biotin Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Biotinylated protein.
Conjugates	Unconjugated
Calculated MW	25kDa/55kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C/1 year

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of Biotinylated Goat IgG, using Biotin Antibody.

Published Papers

Wang Shengqi;Li Jing;Xu Shang;Wang Neng;Pan Bo;Yang Bowen;Zheng Yifeng;Zhang Juping;Peng Fu;Peng Cheng;Wang Zhiyu et al., Baohuoside I chemosensitises breast cancer to paclitaxel by suppressing extracellular vesicle/CXCL1 signal released from apoptotic cells, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.