

GM130 Rabbit mAb

Catalog No: #56019

Package Size: #56019-1 50ul #56019-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GM130 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat Monkey Cow Dog
Specificity	GM130 Antibody detects endogenous levels of total GM130
Immunogen Description	A synthesized peptide derived from human GM130
Conjugates	Unconjugated
Other Names	GM130; Gm130 autoantigen; GOLGA 2; Golga2; Golgi autoantigen; Golgin 95;
Accession No.	Uniprot:Q08379
Calculated MW	130kDa
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

IHC:1:50~1:200


ICC/IF:1:50~1:200

IP:1:50


Images




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Immunofluorescent analysis using the Antibody at 1:500 dilution.(HeLa)



Immunohistochemical analysis of paraffin-embedded human cervix carcinoma, using GM130 Antibody.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.