

ErbB2 (HER2) Rabbit mAb

Catalog No: #56035

Package Size: #56035-1 50ul #56035-2 100ul

Description

Product Name	ErbB2 (HER2) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	ErbB2 (HER2) Antibody detects endogenous levels of total ErbB (HER2)
Immunogen Description	A synthesized peptide derived from human ErbB2 (HER2)
Conjugates	Unconjugated
Other Names	C-erbB-2; ErbB2; HER2, MLN 19; NEU; NEU proto-oncogene; NGL; Receptor protein-tyrosine kinase erbB-2;
Accession No.	Uniprot:P04626
Calculated MW	185kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.