

SOX2 Rabbit mAb

Catalog No: #56043

Package Size: #56043-1 50ul #56043-2 100ul


Description

Product Name	SOX2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse
Specificity	SOX2 Antibody detects endogenous levels of total SOX2
Immunogen Description	A synthesized peptide derived from human SOX2
Conjugates	Unconjugated
Other Names	SOX2; ANOP3; MCOPS3; MGC2413; lcc; SRY related HMG box 2;
Accession No.	Uniprot:P48431
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:1000~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200


Images



Western blot analysis of SOX2 expression in (1) NCCIT cell lysate; (2) F9 cell lysate.



Immunohistochemical analysis of paraffin-embedded human cervix carcinoma, using SOX2 Antibody.



Immunofluorescent analysis of NCCIT cells, using SOX2 Antibody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.