

GATA2/3 Rabbit mAb

Catalog No: #56082

Package Size: #56082-1 50ul #56082-2 100ul

Description

Product Name	GATA2/3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	GATA2/3 Antibody detects endogenous levels of total GATA2/3
Immunogen Description	A synthesized peptide derived from human GATA2/3
Conjugates	Unconjugated
Other Names	DCML;GATA binding protein 2;GATA binding protein 3;HDR;HDRS;IMD21;MONOMAC;NFE1B
Accession No.	Uniprot:P23769/P23771
Calculated MW	51kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

IHC:1:50~1:200


ICC/IF:1:50~1:200

IP:1:50

Images



Western blot analysis of GATA2/3 expression in K562 cell lysate.




Immunohistochemical analysis of paraffin-embedded mouse kidney, using GATA2/3 Antibody.

GATA2/3 Antibody (green)



Immunofluorescent analysis of K562 cells, using GATA2/3 Antibody .



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.