

Cyclin A1/A2 Rabbit mAb

Catalog No: #56107

Package Size: #56107-1 50ul #56107-2 100ul

Description

Product Name	Cyclin A1/A2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Cyclin A1/A2 Antibody detects endogenous levels of Cyclin A1/A2
Immunogen Description	A synthesized peptide derived from human Cyclin A1/A2
Conjugates	Unconjugated
Other Names	CCN1; CCNA; CCNA1; CCNA2; CT146; Cyclin-A; Cyclin-A1; Cyclin-A2;
Accession No.	Uniprot:P20248/P78396
Calculated MW	52,49kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000

IHC:1:100~1:500

IP:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.