

Cdc25C Rabbit mAb

Catalog No: #56114

Package Size: #56114-1 50ul #56114-2 100ul

Description

Product Name	Cdc25C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	Cdc25C Antibody detects endogenous levels of total Cdc25C
Immunogen Description	A synthesized peptide derived from human Cdc25C
Conjugates	Unconjugated
Other Names	CDC25; CDC25C; Cell division cycle 25 homolog C; MPIP3; PPP1R60;
Accession No.	Uniprot:P30307
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:10~1:500

ICC/IF:1:50~1:200

IP:1:30

FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.