

TGF beta 1 Rabbit mAb

Catalog No: #56126

Package Size: #56126-1 50ul #56126-2 100ul


Description

Product Name	TGF beta 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TGF beta 1 Antibody detects endogenous levels of total TGF beta 1
Immunogen Description	A synthesized peptide derived from human TGF beta 1
Conjugates	Unconjugated
Other Names	CED; diaphyseal dysplasia 1, progressive; DPD1; TGF-beta 1; TGFB1;
Accession No.	Uniprot:P01137
Calculated MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

Images



Western blot analysis of TGF beta 1 expression in A549 lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.