

Cathepsin L/V/K/H Rabbit mAb

Catalog No: #56147

Package Size: #56147-1 50ul #56147-2 100ul

Description

Product Name	Cathepsin L/V/K/H Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1633
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Cathepsin L/V/K/H Antibody detects endogenous levels of total Cathepsin L/V/K/H
Immunogen Description	A synthesized peptide derived from human Cathepsin L/V/K/H
Conjugates	Unconjugated
Other Names	Cathepsin; CTSL; CTSV; CTSK; CTSH;
Accession No.	Uniprot:P07711/O60911/P43235/P09668
Calculated MW	38kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

IHC:1:50~1:200

ICC/IF:1:50~1:200


FC:1:50

Images



Western blot analysis of Cathepsin L/V/K/H expression in (1) HepG2 cell lysate; (2) NIH/3T3 cell lysate; (3) PC-12 cell lysate;

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.