

Clathrin heavy chain Rabbit mAb

Catalog No: #56165

Package Size: #56165-1 50ul #56165-2 100ul

Description


Product Name	Clathrin heavy chain Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Clathrin heavy chain Antibody detects endogenous levels of total Clathrin heavy chain
Immunogen Description	A synthesized peptide derived from human Clathrin heavy chain
Conjugates	Unconjugated
Other Names	CHC17; Clathrin heavy chain; Clathrin heavy chain like 1; CLH17; CLH22; CLTC; CLTCL; CLTD;
Accession No.	Uniprot:Q00610/P53675
Calculated MW	191kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:1000~1:5000

IHC:1:50~1:200

Images



Western blot analysis of Clathrin heavy chain expression in Hela cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.