

CHD4 Rabbit mAb

Catalog No: #56169

Package Size: #56169-1 50ul #56169-2 100ul

Description

Product Name	CHD4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	CHD4 Antibody detects endogenous levels of total CHD4
Immunogen Description	A synthesized peptide derived from human CHD4
Conjugates	Unconjugated
Other Names	ATP dependent helicase CHD4; CHD4; Chromodomain helicase DNA binding protein 4; Mi 2 autoantigen 218 kDa protein; Mi 2b; Mi2 beta;
Accession No.	Uniprot:Q14839
Calculated MW	260kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000ICC/IF:1:50~1:200FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.