

CD276 Rabbit mAb

Catalog No: #56177

Package Size: #56177-1 50ul #56177-2 100ul


Description

Product Name	CD276 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CD276 Antibody detects endogenous levels of total CD276
Immunogen Description	A synthesized peptide derived from human CD276
Conjugates	Unconjugated
Other Names	4lg B7 H3; B7 homolog 3; B7H3; B7RP-2; CD276 antigen; CD276 molecule;
Accession No.	Uniprot:Q5ZPR3
Calculated MW	Predicted band size: 57 kDa
SDS-PAGE MW	Observed band size: 100 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes: CD276 Rabbit mAb at 1/1k dilution

Lane 1 : 293T whole cell lysates

Lysates/proteins at 20 μ g per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 57 kDa

Observed band size: 100 kDa

Exposure time: 6 seconds



All lanes :CD276 Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HAP1 cell lysate

Lane 2 : CD276 knockdown HAP1 cell lysate

Lysates/proteins at 20 µg per lane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.