

Atg4B Rabbit mAb

Catalog No: #56184

Package Size: #56184-1 50ul #56184-2 100ul

Description


Product Name	Atg4B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Atg4B Antibody detects endogenous levels of total Atg4B
Immunogen Description	A synthesized peptide derived from human Atg4B
Conjugates	Unconjugated
Other Names	Cysteine protease ATG4B; Autophagin-1; Autophagy-related cysteine endopeptidase 1; Autophagy-related protein 4 homolog B; hAPG4B;
Accession No.	Uniprot:Q9Y4P1
Calculated MW	44kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

FC:1:100

Images



Western blot analysis of Atg4B expression in Ramos cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.