

DDR2 Rabbit mAb

Catalog No: #56221

Package Size: #56221-1 50ul #56221-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DDR2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	DDR2 Antibody detects endogenous levels of total DDR2
Immunogen Description	A synthesized peptide derived from human DDR2
Conjugates	Unconjugated
Other Names	CD167b; NTRKR3; Receptor related 3; TKT;
Accession No.	Uniprot:Q16832
Calculated MW	97kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of DDR2 expression in (1) Jurkat cell lysate; (2) NIH/3T3 cell lysate.

Published Papers

el at., Overcoming Matrix Barrier for Enhanced Immune Infiltration Using Metal Organic Framework Coated with siRNA, , (2024)

PMID:

Cheng Zeng;Xiaojing Chen;Mingxi Lin;Yizi Jin;Qing Guo;Teng Zhou;Xingang Wang;Yiping Li;Xinghui Wang;Yongming Han;Ling Du;Qianyun

Tang;Peifeng Liu;Jian Zhang et al., Overcoming matrix barriers for enhanced immune infiltration using siRNA-coated metal-organic frameworks., , (2025)

[PMID:40054648](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.