

AMBP Rabbit mAb

Catalog No: #56224

Package Size: #56224-1 50ul #56224-2 100ul

Description

Product Name	AMBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	AMBP Antibody detects endogenous levels of total AMBP
Immunogen Description	A synthesized peptide derived from human AMBP
Conjugates	Unconjugated
Other Names	A1M; alpha 1 microglobulin bikunin; Alpha 1 microglycoprotein; AMBP; Bikunin; EDC1; HI30;
Accession No.	Uniprot:P02760
Calculated MW	39kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.