

Dynamin 1 Rabbit mAb

Catalog No: #56249

Package Size: #56249-1 50ul #56249-2 100ul

Description

Product Name	Dynamin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Dynamin 1 Antibody detects endogenous levels of total Dynamin 1
Immunogen Description	A synthesized peptide derived from human Dynamin 1
Conjugates	Unconjugated
Other Names	B dynamin; D100; DNM 1; DNM1; Dynamin; Dynamin1;
Accession No.	Uniprot:Q05193
Calculated MW	97kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:50~1:200

ICC/IF:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.