

Alkaline phosphatase Rabbit mAb

Catalog No: #56377

Package Size: #56377-1 50ul #56377-2 100ul

Description

Product Name	Alkaline phosphatase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	intestinal alkaline phosphatase Antibody detects endogenous levels of total intestinal alkaline phosphatase
Immunogen Description	A synthesized peptide derived from human intestinal alkaline phosphatase
Conjugates	Unconjugated
Other Names	Alpi; Glycerophosphatase; IAP;
Accession No.	Uniprot:P09923
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.