

## ABHD5 Rabbit mAb

Catalog No: #56401

Package Size: #56401-1 50ul #56401-2 100ul

## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | ABHD5 Rabbit mAb                                                                                   |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC FC                                                                                          |
| Species Reactivity    | Human                                                                                              |
| Specificity           | ABHD5 Antibody detects endogenous levels of total ABHD5                                            |
| Immunogen Description | A synthesized peptide derived from human ABHD5                                                     |
| Conjugates            | Unconjugated                                                                                       |
| Other Names           | ABHD5; Abhydrolase domain containing 5; CDS; CGI58; IECN2; NCIE2;                                  |
| Accession No.         | Uniprot:Q8WTS1                                                                                     |
| Calculated MW         | 45kDa                                                                                              |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                   |


## Application Details

WB:1:500~1:2000


IHC:1:100~1:500

FC:1:100

## Images



Western blot analysis of ABHD5 expression in HepG2 cell lysate.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.