

eIF4A1 Rabbit mAb

Catalog No: #56446

Package Size: #56446-1 50ul #56446-2 100ul

Description

Product Name	eIF4A1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	eIF4A1 Antibody detects endogenous levels of total eIF4A1
Immunogen Description	A synthesized peptide derived from human eIF4A1
Conjugates	Unconjugated
Other Names	DDX2; DDX2A; eIF 4A I; EIF4A; eIF4A I; eIF4A-I; EIF4A1;
Accession No.	Uniprot:P60842
Calculated MW	46kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.