

CLASP1 Rabbit mAb

Catalog No: #56478

Package Size: #56478-1 50ul #56478-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLASP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Mouse, Rat, Human
Specificity	CLASP1 Antibody detects endogenous levels of total CLASP1
Immunogen Description	A synthesized peptide derived from human CLASP1
Conjugates	Unconjugated
Other Names	clasp1; hOrbit1; MAST1;
Accession No.	Uniprot:Q7Z460
Calculated MW	169kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000

IHC:1:50~1:200

ICC/IF:1:50~1:200

FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.