

# Aminopeptidase A / CD249 Rabbit mAb

Catalog No: #56486

Package Size: #56486-1 50ul #56486-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | Aminopeptidase A / CD249 Rabbit mAb                                                                |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC                                                                                             |
| Species Reactivity    | Human                                                                                              |
| Specificity           | Aminopeptidase A / CD249 Antibody detects endogenous levels of total Aminopeptidase A / CD249      |
| Immunogen Description | A synthesized peptide derived from human Aminopeptidase A / CD249                                  |
| Conjugates            | Unconjugated                                                                                       |
| Other Names           | Aminopeptidase A; APA; bp1; CD249; EAP; Enpep; Gp160; Ly51;                                        |
| Accession No.         | Uniprot:Q07075                                                                                     |
| Calculated MW         | 160kDa                                                                                             |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                   |

## Application Details


WB:1:1000~1:5000

IHC:1:50~1:200

## Images



Western blot analysis of Aminopeptidase A / CD249 expression in human fetal lung lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.