

GDF3 Rabbit mAb

Catalog No: #56498

Package Size: #56498-1 50ul #56498-2 100ul


Description

Product Name	GDF3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	GDF3 Antibody detects endogenous levels of total GDF3
Immunogen Description	A synthesized peptide derived from human GDF3
Conjugates	Unconjugated
Other Names	ecat9; GDF3; KFS3; MCOP7; Vgr2;
Accession No.	Uniprot:Q9NR23
Calculated MW	16/41kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

Images



Western blot analysis of GDF3 expression in 293T cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.