

AP2S1 Rabbit mAb

Catalog No: #56525

Package Size: #56525-1 50ul #56525-2 100ul

Description

Product Name	AP2S1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	AP2S1 Antibody detects endogenous levels of total AP2S1
Immunogen Description	A synthesized peptide derived from human AP2S1
Conjugates	Unconjugated
Other Names	AP17; AP17 delta; Ap2s1; CLAPS2; Sigma2 adaptin;
Accession No.	Uniprot:P53680
Calculated MW	17 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IP:1:50


FC:1:50

Images



Western blot analysis of AP2S1 expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.