

COPS3 Rabbit mAb

Catalog No: #56528

Package Size: #56528-1 50ul #56528-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	COPS3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	COPS3 Antibody detects endogenous levels of total COPS3
Immunogen Description	A synthesized peptide derived from human COPS3
Conjugates	Unconjugated
Other Names	COP9 subunit 3; cops3; CSN3; JAB 1 containing signalosome subunit 3; SGN3; Signalosome subunit 3;
Accession No.	Uniprot:Q9UNS2
Calculated MW	48 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000

ICC/IF:1:50~1:200


IP:1:50

FC:1:50

Images



Western blot analysis of COPS3 expression in HeLa cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.