

## CD97 Rabbit mAb

Catalog No: #56568

Package Size: #56568-1 50ul #56568-2 100ul

## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | CD97 Rabbit mAb                                                                                    |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC FC                                                                                          |
| Species Reactivity    | Human Rat                                                                                          |
| Specificity           | CD97 Antibody detects endogenous levels of total CD97                                              |
| Immunogen Description | A synthesized peptide derived from human CD97                                                      |
| Conjugates            | Unconjugated                                                                                       |
| Other Names           | ADGRE5; CD97; TM7LN1;                                                                              |
| Accession No.         | Uniprot:P48960                                                                                     |
| Calculated MW         | 55/28kDa                                                                                           |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                   |


## Application Details

WB:1:500~1:2000


IHC:1:50~1:200

FC:1:50

## Images



Western blot analysis of CD97 expression in human tonsil cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.