

Desmoplakin Rabbit mAb

Catalog No: #56575

Package Size: #56575-1 50ul #56575-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Desmoplakin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Rat
Specificity	Desmoplakin Antibody detects endogenous levels of total Desmoplakin
Immunogen Description	A synthesized peptide derived from human Desmoplakin
Conjugates	Unconjugated
Other Names	Desmoplakin I; Desmoplakin II; DPI; DPII; DSP; KPPS2; PPKS2;
Accession No.	Uniprot:P15924
Calculated MW	270kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200


FC:1:50

Images



Western blot analysis of Desmoplakin expression in A431 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.