

ATP7b Rabbit mAb

Catalog No: #56752

Package Size: #56752-1 50ul #56752-2 100ul

Description

Product Name	ATP7b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	ATP7b Antibody detects endogenous levels of total ATP7b
Immunogen Description	A synthesized peptide derived from human ATP7b
Conjugates	Unconjugated
Other Names	ATP7B; PWD; WC1; WD; WND; WND/140 kDa;
Accession No.	Uniprot:P35670
Calculated MW	157kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.