

ARF5 Rabbit mAb

Catalog No: #56812

Package Size: #56812-1 50ul #56812-2 100ul

Description

Product Name	ARF5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	ARF5 Antibody detects endogenous levels of total ARF5
Immunogen Description	A synthesized peptide derived from human ARF5
Conjugates	Unconjugated
Other Names	ADP ribosylation factor 5; Arf5;
Accession No.	Uniprot:P84085
Calculated MW	18kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.