

RanBP9 Rabbit mAb

Catalog No: #56840

Package Size: #56840-1 50ul #56840-2 100ul

Description

Product Name	RanBP9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	RanBP9 Antibody detects endogenous levels of total RanBP9
Immunogen Description	A synthesized peptide derived from human RanBP9
Conjugates	Unconjugated
Other Names	BPM90; BPML; IBAP 1; Imp 9; Importin 9; Ran BP9; Ran-binding protein 9; RanBP7; RanBP9; RanBPM;
Accession No.	Uniprot:Q96S59
Calculated MW	91kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.