

CHMP2B Rabbit mAb

Catalog No: #56844

Package Size: #56844-1 50ul #56844-2 100ul

Description

Product Name	CHMP2B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	CHMP2B Antibody detects endogenous levels of total CHMP2B
Immunogen Description	A synthesized peptide derived from human CHMP2B
Conjugates	Unconjugated
Other Names	ALS17; CHMP2b; DMT1; hVps2-2; VPS2B;
Accession No.	Uniprot:Q9UQN3
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.