

EHMT2 Rabbit mAb

Catalog No: #56905

Package Size: #56905-1 50ul #56905-2 100ul

Description

Product Name	EHMT2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	EHMT2 Antibody detects endogenous levels of total EHMT2
Immunogen Description	A synthesized peptide derived from EHMT2
Conjugates	Unconjugated
Other Names	Bat8; Ehmt2; Euchromatic histone lysine methyltransferase 2; G9a protein; GAT8; Histone H3 K9 methyltransferase 3; NG36;
Accession No.	Uniprot:Q96KQ7
Calculated MW	160,170kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.