

AP2S1 Rabbit mAb

Catalog No: #56927

Package Size: #56927-1 50ul #56927-2 100ul

Description

Product Name	AP2S1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	AP2S1 Antibody detects endogenous levels of total AP2S1
Immunogen Description	A synthesized peptide derived from AP2S1
Conjugates	Unconjugated
Other Names	AP17; Ap2s1; CLAPS2; Sigma adaptin 3b; Sigma2 adaptin;
Accession No.	Uniprot:P53680
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Note: This product is for in vitro research use only and is not intended for use in humans or animals.